
































Matrix Back: Order Assistant



Process Stages:

- 1) Measure Patient, C, or use Back Template Set
- 2) Choose Back size and type
- 3) Specify mounting hardware once wheelchair width and back upright diameter are known; see section 4A.

Size	C	width	Code (BC-x) Back Complete	*	** Basic Back	Code (BCTB-x) Back Complete 2 Balls	*	** Basic Back Two Balls
1	16"	25"	 <p>88 Clamps 116 4-Balls 19 Green 4-Balls 28 Black Balls 14 Half 4-Balls</p> <p>BC-1</p>			 <p>88 Clamps 84 4-Balls 17 Green 4-Balls 28 Black Balls 14 Half 4-Balls 24 2-Balls</p> <p>BCTB-1</p>		
2	17"	27"	 <p>106 Clamps 150 4-Balls 20 Green 4-Balls 32 Black Balls 16 Half 4-Balls</p> <p>BC-2</p>			 <p>106 Clamps 122 4-Balls 18 Green 4-Balls 32 Black Balls 16 Half 4-Balls 30 2-Balls</p> <p>BCTB-2</p>		
3	18"	29"	 <p>127 Clamps 182 4-Balls 22 Green 4-Balls 33 Black Balls 16 Half 4-Balls</p> <p>BC-3</p>			 <p>127 Clamps 148 4-Balls 20 Green 4-Balls 33 Black Balls 16 Half 4-Balls 36 2-Balls</p> <p>BCTB-3</p>		
4	19"	21"	 <p>142 Clamps 216 4-Balls 23 Green 4-Balls 38 Black Balls 19 Half 4-Balls</p> <p>BC-4</p>			 <p>142 Clamps 182 4-Balls 21 Green 4-Balls 38 Black Balls 19 Half 4-Balls 36 2-Balls</p> <p>BCTB-4</p>		
5	20"	33"	 <p>168 Clamps 258 4-Balls 26 Green 4-Balls 42 Black Balls 14 Half 4-Balls</p> <p>BC-5</p>			 <p>168 Clamps 216 4-Balls 26 Green 4-Balls 42 Black Balls 14 Half 4-Balls 42 2-Balls</p> <p>BCTB-5</p>		

* Complete with all mounting hardware and one black cover. Please complete hardware order form.
 ** Matrix Back Sheet and T-Track ONLY. No Cover, hinges, arms or central assembly hardware.
 C = Accurately measure from ASIS to top of shoulder on patient or use Back Template Set